


Calibration in the CDR FOODLAB® line analyzers: Design and Flexibility

All the chemical analysis systems in **CDR FOODLAB® line** are engineered to ensure precision, reliability, and reproducibility through factory calibration and built-in self-check procedures that maintain long-term analytical consistency. This article explores the calibration process, why it is typically not required, and how the system's optional calibration features provide advanced users with the flexibility to align results with alternative methods or historical data sets.

CDR FOODLAB® Calibration: Simple, Reliable, Flexible

Precision Out of the Box



Instruments and reagents are factory-calibrated, saving time and eliminating the need for routine manual setup.

Self-Checks for Lasting Accuracy



At every startup, the system automatically verifies its calibration integrity, ensuring consistent, reliable results.

Align Data Your Way



Optionally fine-tune the system to match historical data or other analytical methods for seamless workflow integration.

1. Factory Calibration: Instant Precision

CDR analytical systems are factory pre-calibrated, ensuring that both reagents and instruments are optimized for accuracy and consistency before reaching the user. The reagents are produced by CDR's chemical laboratories to provide analytical results aligned with official reference methods. Each new batch is prepared and calibrated to maintain full consistency with previous lots, eliminating the need for any further adjustment. Internal quality control monitors the accuracy of every production step, from reagent preparation to the filling of each cuvette, ensuring reproducibility and long-term reliability. Likewise, during production, each instrument undergoes calibration of the optical reading unit, performed according to a standardized protocol. This process guarantees that all CDR analysis systems worldwide deliver the same analytical response and operate in a perfectly consistent way.

Calibration data are stored directly within the analyzer, allowing the system to be immediately ready for use without the need for routine manual calibration. This integration between instrument and reagents ensures high analytical performance, saves time and costs, and minimizes the potential for operator error.

2. Automated Self-Check: Ensuring Ongoing Reliability

To preserve the integrity of the calibration performed during production, CDR systems are equipped with an automatic verification mechanism for analytical reliability that is activated at each startup. This internal diagnostic process verifies that the stored calibration parameters are intact and that the system is operating within its specified performance range. If any discrepancies are detected, the system alerts the user, ensuring that



analisi e sviluppo
sistemi cibernetici

CDR S.r.l. Via degli Artigiani, 3-6
50055 Ginestra F.na | Firenze
Tel. 055.871431 | Fax 055.8714322
www.cdr-mediated.it

Reg. Imprese Fi e C.F. n. 04170450482
P. IVA IT 04170450482 - R.E.A. FI 423250
Cap. Soc. € 60.000 - i.v.



measurements are not compromised. This self-check feature reinforces the system's reliability, making it suitable for environments where consistent results are critical, such as quality control, research, or regulatory compliance. This design philosophy prioritizes ease of use while maintaining the precision required for professional applications.

3. Optional Fine Tuning: Flexibility for Specialized Needs

Despite calibration performed during production and automatic verification ensuring high precision, CDR systems also provide calibration functions for users with specific requirements. These functions allow the system to be fine-tuned to align results with those obtained using alternative analytical methods or to ensure consistency with historical data series, providing a high degree of customization. For example, laboratories transitioning from older analytical systems to those of the CDR FOODLAB[®] line may need to ensure that new measurements are consistent with existing data archives. By performing an optical calibration, users can adjust CDR systems to produce results that align with those generated by official methods, facilitating integration into established workflows. Similarly, in cases where a different reference method is used (e.g., a specific chromatographic or spectroscopic technique), the system can be calibrated to harmonize its output with these methods, ensuring comparability across platforms. This optional calibration is particularly valuable in research and industrial settings where data continuity or cross-method validation is essential. The process is user-friendly, guided by intuitive software interfaces, and does not require advanced technical expertise.

4. Benefits of the CDR Calibration Approach

The approach of CDR analysis systems to calibration offers several key benefits:

- **Simplicity and Efficiency:** Factory calibration and automated self-checks eliminate the need for routine manual calibration, reducing setup time and operational complexity.
- **Reliability:** The self-check mechanism ensures that the system remains within calibration specifications, making each measurement reliable.
- **Flexibility:** Optional calibration allows users to tailor the system to specific needs, such as aligning with alternative methods or historical data, without compromising ease of use.
- **Versatility:** The system supports a wide range of applications, from food and beverage analysis to environmental monitoring, with the ability to adapt to diverse laboratory requirements.
- **Cost-Effectiveness:** By minimizing the need for external calibration tools or frequent recalibration, the system reduces operational costs while maintaining high performance.

5. Fine-tuning for data consistency

The possibility of fine-tuning CDR analyzers is a feature which is particularly useful in scenarios where data consistency across different systems or time periods is critical. For instance:

- **Quality Control in Food Production:** A food manufacturer may need to align the CDR system's results with those of an older instrument used in previous production cycles to ensure compliance with historical quality standards.
- **Ensuring Data Continuity in Quality Control:** In beverage production, comparing new analytical data with historical results is essential to understand how processes evolve over time. By calibrating the CDR analyzer to match previously used systems, manufacturers can maintain data consistency, enabling accurate long-term trend analysis and continuous improvement of quality control.



analisi e sviluppo
sistemi cibernetici

CDR S.r.l. Via degli Artigiani, 3-6
50055 Ginestra F.na | Firenze
Tel. 055.871431 | Fax 055.8714322
www.cdr-mediated.it

Reg. Imprese Fi e C.F. n. 04170450482
P. IVA IT 04170450482 - R.E.A. FI 423250
Cap. Soc. € 60.000 - i.v.



- **Method Integration:** Laboratories that adopt the CDR analyzer alongside other analytical techniques (e.g., HPLC or titration) can use calibration to ensure that results are comparable, enabling effective integration of analytical workflows.

Conclusion

CDR analytical systems are designed to deliver accurate and reliable results without any user intervention, thanks to production-calibrated reagents and instruments combined with automatic verification of analytical reliability at startup. This streamlined approach makes the systems ideal for users seeking simplicity and consistency. For those with specialized needs, the optional calibration feature provides the flexibility to align results with alternative methods or historical data, ensuring compatibility with existing workflows and data sets. By combining ease of use, reliability, and adaptability, CDR systems stand out as a versatile solution for a wide range of analytical applications.

Summary Table of Key Concepts

Topic	Summary
Factory Calibration	Instruments and reagents are pre-calibrated during production to ensure accuracy, reproducibility and global consistency. No routine calibration is required.
Reagent Calibration	Each reagent batch is produced and validated by CDR's chemical laboratories to match reference methods and maintain continuity with previous batches.
Instrument Calibration	The optical units are calibrated during production using standardized procedures. Calibration data are stored in the analyzer for immediate use.
Automatic verification of analytical reliability	At each startup, the system checks the stored calibration parameters and performance conditions, alerting the user to any deviations in order to ensure consistent reliability.
Optional Fine Tuning	Users can align system outputs with alternative analytical methods or historical datasets. Useful for method transition, research or long-term data continuity.
Use Cases for Fine Tuning	Alignment with previous instruments, maintenance of data continuity in food and beverage production, and support for method integration with techniques such as HPLC or titration.
Benefits	Simplified workflows, no routine calibration, high reliability, customizable alignment with reference methods and historical data, a versatile range of applications, and reduced operating costs.
Overall Conclusion	CDR systems deliver precise and consistent results without any intervention, while also offering flexibility when specific calibration alignment is required.



analisi e sviluppo
sistemi cibernetici

CDR S.r.l. Via degli Artigiani, 3-6
50055 Ginestra F.na | Firenze
Tel. 055.871431 | Fax 055.8714322
www.cdr-mediated.it

Reg. Imprese Fi e C.F. n. 04170450482
P. IVA IT 04170450482 - R.E.A. FI 423250
Cap. Soc. € 60.000 - i.v.

