



CDR FoodLab[®], the system for the quality control of Oils&Fats

SIMONE BELLASSAI CDR FoodLab[®] Division Manager



CDR, a "workshop" of ideas and continuous innovation

CDR conducts its business in heterogeneous sectors



TELEMATIC SYSTEMS Automatic Toll collection terminals

MEDICAL DIAGNOSTICS Hematology and hemostasis systems

FOOD DIAGNOSTICS Analysis systems for food and beverage SENSORING SYSTEMS Sensors and probes for QC

SECTORS









- Milk and Dairy Products
- Egg Products
- Tomato/Vegetable puree
- Bakery products



- Wine
- Beer and Water
- Cider
- Kombucha
- CDR WineLab[®]
- CDR BeerLab[®]
- CDR CiderLab
- CDR KombuchaLab



- All kinds of vegetable oils
- All kinds of animal fats
- Nuts/hard-shelled fruits



- CDR OxiTester
- CDR PalmOilTester





Canada USA Mexico Guatemala Jamaica Trinidad e Tobago

> Colombia Ecuador Suriname Perù Chile Argentina Brazil Paraguay

> > Uruguay

More than 100 countries

Estonia Hungary Island Latvia Romania Finland Lithuania Bulgaria Norway Poland Serbia Sweden Cezh Republic Slovenia Denmark Ukraine Croatia Russsia Slovakia Albania Austria Switzerland Ireland Belgium France United Kingdom Netherland Germany Monaco Spain Portugal Greece Morocco Algeria Tunisie Libya Egypt Sudan Seychelles Uganda Kenya **Ivory Coast** Tanzania Togo Burundi Ghana Ruanda Nigeria Botswana Cameroon South Africa Gabon Congo

Italy

Turkey Azerbaijan Cyprus Syria Irag Kuwait Iran Jordan Lebanon Israel Palestine Saudi Arabia Emirates Baharain

China Mongolia Japan South Korea Philippines Mauritius Myanmar Vietnam Thailand Taiwan Bangladesh India Pakistan Indonesia Malaysia Singapore Sri Lanka

Australia New Zeland New Caledonia Papua New Guinea



CDR FoodLab[®] is a range of chemical analysis systems, easy to use and versatile, that allows to determine a wide panel of parameters on food and beverage



CDR FoodLab[®] is composed by thermostated analyzers based on photometric technology that uses LED emitters at fixed wavelengths



Analysis kits are supplied in package of 10 tests, containing ready to use pre-filled reagents, designed and produced by CDR



Which are the most important features and advantages of CDR systems compared to the traditional methods?



Easy to use





Possibility to use the analyzer at line in the processing plant

Maintenance free

Calibration free

Everyone can use CDR FoodLab[®] line



INSTRUMENT's FEATURES

- ✓ No maintenance
- ✓ No calibration!!!
- ✓ 16 analyses in about 10 mins (analyses by session)
- ✓ Multitasking mode (different analyses at the same time)
- ✓ 3 years warranty
- Touch screen where procedures are explained step by step

Reading cells and incubation positions are thermostated at 37°C





INSTRUMENT's FEATURES

The reading cell



Standard photometry

The range in absorbance is 0,000 - 6,000

The range in absorbance is 0,000 - 2,000/3,000

The range in absorbance is 2-3 times more than normal standard photometry



CDR FoodLab[®] system



CDR FoodLab®

It can do analysis on all the matrices available

- With printer
- 16 analysis of the same type at the same time
- MULTITASKING mode that means different a analysis at the same time



CDR FoodLab[®] Jr.

It can be configured with FFA, PV, AnV, Soaps and Iodine Value

- Without printer
- 3 analysis of the customer choice
- No MULTITASKING



REAGENTS' FEATURES

- ✓ Pre-filled cuvettes
- ✓ Ready to use reagents
- ✓ Microquantity of sample
- ✓ You don't need any calibration
- ✓ Long shelf-life of the reagents





THE READING CELL



Standard photometry

The detectable range of ABSORBANCE is 0,000 – 6,000

The detectable range of ABSORBANCE is 0,000 – 2,000/3,000

The range of ABSORBANCE of CDR FoodLab[®] is 2-3 times wide as it is wide the one of a standard spectrophotometer



Comparison with standard photometry



CDR reagents CALIBRATED AND READY TO USE

- Thermostated instrument
- Short reaction times thanks to thermostated element
- High preformance reading cell in terms of calibration and long-life
- Maintenance not required
- Skilled technicians not required
- Wider ABSORBANCE range (ABS 0,000 6,000)

Standard photometry

- Enzymatic kits to be prepared and calibrated
- Instrument usually not thermostated
- Longer reaction times
- Short-life reading cell, in case of frequent use
- Maintenance required
- Skilled technicians required
- Shorter ABS range (ABS 0,000 2,000)



ADVANTAGES of CDR FoodLab[®]

- Easy to use
- Strong system
- No calibration or check is required
- Everyone can use the system, you don't need any technical background
- Wider Absorbance range (ABS 0,000 6,000)
- Possibility to get real time results



CDR FoodLab[®] system



With CDR FoodLab[®]you can carry out easily and quickly on different matrices analysing **liquid** and **solid samples**



CDR FoodLab[®] system

Bakery products

- Free Fatty Acids (FFA)
- Peroxide Value
- p-Anisidine Value
- Lactose
- Alcohol on final product
- Alcohol by Volume in the antibacterial solution



Egg Products

- Lactic Acid
- D-3-hydroxybutyric acid
- Colour (Beta-Carotene)
- Cholesterol



Milk and Dairy products

- Lactose
- L-Lactic acid
- Milk Urea Nitrogen (MUN)
- Alkaline phosphatase (ALP)
- Ammonia
- Chloride
- Hydrogen peroxide
- ε-fructosyl-lisine (Furosine)
- Peroxidase

Oils&Fats

- Free Fatty Acids
- Peroxide Value
- Anisidine Value (AnV)
- Soaps Test
- Iodine Value test

Tomate

- Reducing Sugar
- Lactic acid (D+L)
- Chloride



You can carry out analysis of FFA, PV and AnV on **Nuts/Hard shelled fruits**



A focus on the most important analyses in Oils&Fats industry

• Free Fatty Acid (FFA)



• Peroxide Value (PV)



It is the most important parameter related to the quality of the oil/fat

It is the primary oxidation state of the oil/fat

It is the secondary oxidation state of the oil/fat



Which are the traditional methods for FFA, PV and AnV?

Free Fatty Acids



Titration

Peroxide Value



Anisidine Value

Chemical reaction detected photometrically



Free Fatty Acids (FFA)

Traditional method (Titration)

- 50-100mL of solvent
- Sodium hydroxide
- Glassware
- Fume hood
- Phenolphthalein

Estimated time: 5 minutes

CDR FoodLab[®] system

- CDR FoodLab[®] analyser
- 1 mL in pre-filled cuvette

Estimated time: 10 seconds

Analyses session of 16 samples at the same time in 10 minutes



Peroxides Values (PV)

Traditional method (Titration)

CDR FoodLab[®] system

- 25mL of mixture glacial acetic acid/chloroforme
- $Na_2S_2O_3$ solution
- KI solution
- Starch solution
- Glassware
- Fume hood

Estimated time: 10 minutes



- CDR FoodLab[®] analyser
- 1 mL in pre-filled cuvette
- 10uL of R2 reagent

Analysis time: 3 minutes

Analyses session of 16 samples at the same time in 10 minutes



Anisidine Value (AnV)

Photometric analysis

- 25mL of isooctane
- Glacial acetic acid
- p-Anisidine in glacial acetic acid
- Photometer (reading at 350nm)
- Glassware
- Fume hood

Estimated time: 15 minutes



CDR FoodLab[®] system

- CDR FoodLab[®] analyser
- 1 mL in pre-filled cuvette

NO carcinogenic reagent!

Analysis time: 1 minute

Kinetic reaction: samples have to be prepared one by one



CDR FoodLab[®] Measuring ranges



• Free Fatty acids (0,01 – 26% oleic acid)

• Peroxides Value $(0,01 - 550 \text{ meqO}_2/\text{Kg})$

• Anisidine Value (0,5 – 100 AnV)



CDR FoodLab[®]

Correlation with the reference method









CDR FoodLab[®] Repeatability

Anisidine Value

	CDR FoodLab® data	Campden BRI reference method data
Replicate number	Anisidine (AnV)	Anisidine (AnV)
1	5.1	5.1
2	5.2	5.1
3	5.2	5.1
4	4.9	5.0
5	5.3	5.2
6	5.3	5.1
7	5	5.2
8	5.1	5.3
9	5.3	5.1
10	5.2	5.1
Mean	5.2	5.1
Standard deviation	0.13	0.07

Table 5 Repeatability for Anisidine Value (using Olive Oil 3)

	CDR FoodLab® data	Campden BRI reference method data
Replicate number	Peroxides (mEq O ₂ /kg)	Peroxides (mEq O ₂ /kg)
1	9.22	9.69
2	10.51	9.13
3	10.56	9.57
4	9.33	9.66
5	9.25	9.57
6	9.38	9.41
7	9.3	9.25
8	9.1	9.92
9	9.14	9.15
10	9.31	9.76
Mean	9.51	9.51
Standard deviation	0.52	0.27

Table 6 Repeatability for Peroxide Value (using Olive Oil 3)

Free Fatty Acid content			
	CDR FoodLab® data	Campden BRI reference method data	
Replicate number	FFA (% oleic acid)	FFA (% oleic acid)	
1	0.36	0.37	
2	0.34	0.35	
3	0.33	0.36	
4	0.34	0.35	
5	0.34	0.35	
6	0.35	0.36	
7	0.34	0.36	
8	0.35	0.35	
9	0.35	0.35	
10	0.34	0.35	
Mean	0.34	0.36	
Standard deviation	0.01	0.01	
Table 7 Repeatability for Free Fatty Acid content (using Olive Oil 3)			



LIQUID SAMPLES



We can analyze crude or refined

We don't have any problem of colour or viscosity

In case of solid fat, it is enough to melt it to collect it and carry out the analysis





Nuts/Hard shelled fruits/ Seeds

Solid Samples



Flours



Pet Food



Analisis of **snacks** with CDR FoodLab[®]

With the **CDR FoodPress** we can easily and quickly extract the oil from the solid sample

After a short centrifugation we take the the sample and we can carry out all the analysis (FFA, PV and AnV)





Analisis on **nuts and hard shelled fruits** with CDR FoodLab[®]

With the **CDR FoodPress** we can easily and quickly extract the oil from the solid sample

After a short centrifugation we take the the sample and we can carry out all the analysis (FFA, PV and AnV)





Analisis of **flour** with CDR FoodLab[®]

The **flour** is treated with a specific extraction solution to extract the fat from the solid sample

In this case CDR FoodLab[®] need the **Fat %** so that at the end the results of FFA, PV and AnV are recalculated automatically by the analyzer





Analisis of **pet food** with CDR FoodLab[®]

Pet Food has to be grinded to get a flour. After that with a specific extraction solution we extract the fat from the solid sample

In this case CDR FoodLab[®] need the **Fat %** so that at the end the results of FFA, PV and AnV are recalculated automatically by the analyzer







Thank you everyone!



Simone Bellassai CDR FOODLAB division manager