



# **CDR WineLab<sup>®</sup>: An easy-to-use device for the real-time monitoring of your fermentations**

**Simone Bellassai**  
**CDR FoodLab<sup>®</sup> division manager**



CDR, a “workshop” of ideas and continuous innovation

## CDR conducts its business in heterogeneous sectors

### SECTORS



#### **TELEMATIC SYSTEMS**

Automatic Toll  
collection terminals



#### **MEDICAL DIAGNOSTICS**

Hematology and hemostasis  
systems



#### **FOOD DIAGNOSTICS**

Analysis systems for  
food and beverage



#### **SENSING SYSTEMS**

Sensors and probes for QC



**Food and Feed**

- Milk and Dairy Products
- Egg Products
- Tomato/Vegetable puree
- Bakery products



➤ **CDR FoodLab®**



**Beverage**

- Wine
- Beer and Water
- Cider
- Kombucha



- **CDR WineLab®**
- **CDR BeerLab®**
- **CDR CiderLab**
- **CDR KombuchaLab**



**Oils and Fats**

- All kind of oils
- All kind of fats
- Nuts/hard-shelled fruits



- **CDR FoodLab®**
- **CDR OxiTester**
- **CDR PalmOilTester**







CDR WineLab<sup>®</sup> is the sytem composed of an analyzer based on photometric technology and pre-filled ready to use reagents



Also available CDR WineLab<sup>®</sup> Junior version.  
Both are thermostated analyzers  
at 37°C that use LED emitters as light  
sources, to specific wavelengths



Reagents are supplied in bag of 10 pre-filled  
vials, ready to use and “calibrated”



Pipette for the sample collection

# ANALYZER FEATURES



- ✓ No maintenance
- ✓ No calibration of the analyzer is required
- ✓ LEDs wavelengths with an Absorbance range 0.000 – 6.000
- ✓ 16 analyses in the same working session
- ✓ MULTITASKING mode (different tests at the same time)
- ✓ 3 years warranty
- ✓ Touch screen on which methods are explained step by step

INCUBATION  
PART

PRINTER



READING PART  
made by **LEDs wavelengths**

The two parts dedicated to  
incubation and reading  
are thermostated at 37°C

# REAGENTS FEATURES



- ✓ Pre-filled reagents
- ✓ NO CALIBRATION REQUIRED
- ✓ Ready to use reagents
- ✓ Long shelf-life (from 6 to 18 months)



# CDR WineLab<sup>®</sup> system



CDR WineLab<sup>®</sup>

You can do all the analyses of the panel

- Printer on board
- You can carry out 16 analyses at the same time of the same parameter
- You can carry out different analyses at the same time (multitasking mode)



CDR WineLab<sup>®</sup> Jr.

It can be tailored with 3 analyses of your choice with the possibility to add the others later

- No printer
- You can carry out 3 analyses at the same time of the same parameter





Analysis on GRAPES



Analysis on GRAPE JUICE



Analysis on WINE



**We can analyze the whole winemaking process,  
from grape to finished wine**



**Total anthocyanins**

**Total Polyphenols Index (TPI)**

**Total acidity**

**pH**

**Yeast Assimilable Nitrogen (YAN)**

**Free sulfur dioxide**

**Tonality and intensity**

**Glycerol**

**Catechins**

**Acetaldehyde**

**Glucose and Fructose**

**L-Lactic acid**

**Total sulfur dioxide**

**Tannins**

**HCl index**

**Gluconic acid**

**L-Malic acid**

**Calcium**

**Copper**

**Galacturonic acid**

**Acetic acid**

**Polymerised anthocyanins**

**Alcohol by volume**

# Winemaking –red wine

Ripeness of grapes



Before starting the pressing phase you can determine ***Technological ripeness*** and ***Phenolic maturity***

Technological  
ripeness



$$\frac{\text{Sugars } (\frac{g}{L})}{\text{Total Acidity } (\frac{g}{L})}$$

Phenolic  
maturity




With CDR WineLab® it is possible to follow the phenolic maturity of grape measuring the anthocyanins

# Winemaking - red wine

Analysis on grape juice



**Grape juice** is a liquid with a high sugar content (about 200g/L) and a high turbidity, formed during the phases of destemming and crushing of grapes



Yeasts will be added  
to the must for starting  
ALCOHOLIC FERMENTATION

For accurate analyses, it is required to **centrifuge** the must.

If must is fermenting, it is required to **degas** it with an **ultrasonic bath**, and then to centrifuge it.



**Centrifuge**



to clarify the grape juice  
(3 minutes)



**Ultrasonic bath**




to degasing the grape juice  
(2-3 minutes)




# Winemaking - red wine


Analysis on grape juice with CDR WineLab®



Glucose



Fructose



Sucrose

## ***Fermentable sugars (reference method)***

CDR WineLab® detects only the two fermentable sugars - glucose and fructose - with three different kits

Fermentable sugars

➔ 230 g/L

Glucose and fructose

➔ 120 g/L  
110 g/L

} Easy dilution kit is used

Fermentable sugars  
(sucrose included)


➔ 250 g/L  
with addition of  
20 g/L of sucrose

In this case, the dilution/inversion kit is included




# Winemaking – red wine

The end of alcoholic fermentation




# Winemaking - red wine


Analysis on grape juice with CDR WineLab®



Tartaric acid



L-Malic acid



Acetic acid


## ***Total acidity (method correlated to the reference one)***

It detects the content of organic acids in must/wine. It is a very important parameter for the wine taste and it can be adjusted by the addiction of tartaric acid to the must.


CDR WineLab® detects the total acidity in the range 0.5 – 10.0 g/L; the result is expressed as tartaric acid  
*CDR WineLab® configuration for France foresees the result as sulfuric acid.*

# Winemaking - red wine


Analysis on grape juice with CDR WineLab®



Acetic acid



L-Malic acid



Sulphur dioxide

## ***Acetic acid (reference method)***

Parameter referred to the presence of bacterial contamination. Its increasing represents a not recoverable organoleptic alteration. Its detection is very important in order to monitor the whole winemaking process.

CDR WineLab® determines acetic acid in the range 0.05 – 1.20 g/L.

## ***L-Malic acid (reference method)***

Fundamental parameter for the malolactic fermentation process. This analysis can be done on the pre-fermented must also, in order to plan the proper fermentation protocol.

CDR WineLab® determines L-malic acid in the range of 0.05 – 5.00 g/L.

## ***Total sulfur dioxide***

It is added before fermentation in order to reduce acetic acid bacteria and the development of indigenous and wild yeasts as much as possible, and to protect the grape juice from oxidation.


# Winemaking - red wine

Analysis on grape juice with CDR WineLab®




## ***Yeast assimilable nitrogen (reference method)***

Parameter that represents the concentration of assimilable nitrogen from yeasts during the fermentation process. It is composed of inorganic nitrogen (ammoniacal nitrogen) and organic nitrogen ( $\alpha$ -amino nitrogen).



Ammonium



Aminoacid



CDR WineLab® allows to detect both nitrogen components separately, in order to optimize the addition of nutrients to the must




Proline and hydroxy-proline acids, widely present in must but not fermentable by yeasts, are not detected.

On the contrary, the old method (titration with formaldehyde) is affected by this interference.

## ***Gluconic acid (reference method) and Glucuronic acids***

Parameter referred to the action of *Botrytis Cynerea* on grapes. The increasing of its concentration in the determines a more difficult protection of the wine as finished product, because gluconic acid is related to an increase of substances that combine with  $\text{SO}_2$



Fundamental parameter for the production of Champagne/Spumante, widely used by our customers located in the region of Champagne: Moët Chandon, Veuve Clicquot and Laboratoire du Champagne




# Winemaking – red wine

## Maceration



The maceration occurs during the alcoholic fermentation and can be analyzed by:

- **Color analysis (ABS 420, 520, 620nm)**
- **Total Anthocyanins**
- **Total Polyphenols Index**



**Centrifuge**




**Ultrasonic bath**



# Winemaking – red wine

## Malolactic fermentation



L-Malic acid test is used for detecting the end of the malolactic fermentation  
CDR WineLab® detection limit is 0.05 g/L

L-Lactic acid test - not L-Malic acid test -  
has to be used for detecting the beginning of the malolactic fermentation!!!

— L-malic acid  
— L-lactic acid

# Winemaking – red wine

Color stabilization and aging



**Barrique and Tonneau**  
(aging in wood)



**Concrete vessels**  
(aging in concrete)



**Steel vessels**  
(MICRO-OXIGENATION)



**Steel vessels**  
(aging in steel)

## Polyphenols monitoring

- **Total anthocyanins**
- **Total Polyphenols Index (TPI)**
- **Total Polyphenols Folin Ciocalteu**
- **Color**



# CDR WineLab® - PARSEC



- Total anthocyanins
- Polymerised anthocyanins
- Total Polyphenols Index (TPI)
- Tannins
- HCl index
- Color
- pH

# Winemaking - white wine

Ripeness of grapes



Before starting the real mashing phase, standard analyses for checking the technological ripeness of grapes are performed

Technological ripeness



$$\frac{\text{Sugars } (\frac{g}{L})}{\text{Total Acidity } (\frac{g}{L})}$$

Most of the white wines are sensible to grey mould, due to the development of *Botrytis Cinerea*. The preventive analysis of **gluconic acid** and **galacturonic acid**, allows the oenologist to better evaluate the health condition of grapes

The quantification of gluconic acid and galacturonic acid are very important parameters in order to evaluate the quality of white grape juice that will become *champagne* or *spumante*





# Winemaking - white wine

## grape juice extraction

All the operative steps before fermentation (on grapes and must) are of huge importance for the evaluation of the product quality. It is necessary to:

- gently press the grapes
- reduce the mechanical pressure that can affect the skins
- increase the pressure in a slow progressive way
- do not apply a temperature over 20°C for extracting the juice
- reduce the turning over of pomace as much as possible
- keep away from air (from oxygen in particular) the extracted must

Checking during extraction phase  
with CDR WineLab®



The must browning can be evaluated  
by **ABS** reading to **420nm**



The extraction of phenolic  
components can be evaluated by  
the **Total Polyphenols Index** test





# Winemaking - white wine

## Fermentation

The evolution of alcoholic fermentation is usually slower in white wines than in red wines, because in this case the temperature is kept lower, in order to better preserve the aromas.

As far as the sugars detection with CDR WineLab® is concerned, this is performed as already explained for red wines.

The malolactic fermentation usually is not done in this case, because L-Malic acid gives acidity and freshness to a white wine.

**Malolactic kit** – allows the determination of both L-Malic and L-Lactic acids

It is used by those winery that produce white wines in which the malolactic fermentation is partially performed

*Example of result:*

1,21 g/L of L-Malic acid


0,45 g/L of L-Lactic acid



# Wine bottling

At the end of the winemaking process for the bottling  $\text{SO}_2$  free and total are two very important parameters

- **Free  $\text{SO}_2$**   
It represents the real protection, from both microbiological and chemical points of view, for wine
- **Total  $\text{SO}_2$**   
Related to the legal limit




**CDR WineLab<sup>®</sup>** for Total and Free  $\text{SO}_2$  uses a reagent not affected neither by wine color nor ascorbic acid


# CDR WineLab® correlation with the reference method




**Acetic acid**



**L-malic acid**



**Fermentable sugars**





## In conclusion

- CDR WineLab<sup>®</sup> system is an easy and fast tool for your wine-making QC
- You can take decisions quickly in a few minutes about the wine making process
- You can realize a complete in-house quality control of the process
- The analyzer can be used by everyone. You don't need any chemical expertise
- You don't need any glassware. With only a small desk you can check the whole production process



**Let's see right now how  
CDR WineLab<sup>®</sup> works!!!**



**Simone Bellassai  
CDR FOODLAB<sup>®</sup> division manager**