

## Acetic Acid Measurement: A Strategic Guide to Wine Quality and Analytical Innovation

Dr. Francesca Bruni, researcher at the CDR Chemical Lab “Francesco Bonicolini”

P2609

### 1. The Role of Volatile Acidity in the Wine Profile

Volatile acidity (VA) is not simply an analytical parameter, but rather an indicator of the microbiological health and stability of wine. Strategically speaking, monitoring this parameter means managing risk: an uncontrolled increase not only compromises the sensory profile, but can lead to the complete loss of the product's commercial value due to acetification.

Chemically, volatile acidity consists primarily of **acetic acid**, which is highly volatile. When its concentration exceeds certain thresholds, the organoleptic impact becomes critical, manifesting itself with the typical "pungent" and acrid aroma that masks the noble aromatic components, a prelude to the wine's degradation into vinegar. The legal and sensory references follow rigorous standards defined by the OIV:

- **Sensory Threshold:** Generally perceptible around **0.5 g/L**, although the structure of the wine may influence its detection.
- **Legal Limits (OIV):** The maximum limit is set at **1.2 g/L** of acetic acid for still white and red wines, which can be raised to **2.1 g/L** for special wines such as straw wines or botrytized wines.
- **Quality Standards:** For excellent productions, the target values are between **0.2 and 0.6 g/L**.

Constant monitoring is the only way to prevent irreparable defects. The adoption of rapid and accurate analytical methods allows for the development of an active prevention strategy, essential for protecting company investments.

### 2. Formation Dynamics: Why and When Acetic Acid Develops

Identifying the source of acetic acid production is essential to distinguish between physiological fermentation dynamics and destructive bacterial contamination, guiding corrective interventions before the damage is irreversible.

The causes of formation are divided into three critical levels:

1. **Metabolism of yeasts and lactic acid bacteria:** During alcoholic fermentation, *Saccharomyces cerevisiae* physiologically produces 0.2–0.4 g/L of acetic acid. In malolactic fermentation (MLF), bacteria such as *Oenococcus oeni* can increase this

value (typically around 0.4 g/L), with risks of peaks above 1.0 g/L in uncontrolled environmental conditions.

2. **Deteriorated grapes and Botrytis:** The health of the raw material is the primary risk factor. While the genus *Gluconobacter* (which prefers glucose and has a low tolerance for acetic acid) predominates on healthy grapes, the presence of ***Acetobacter increases on grapes damaged or affected by Botrytis cinerea***. These bacteria consume the ethanol produced by yeasts to generate acetic acid.
3. **Management deficiencies and oxygen:** Acetic acid is the product of aerobic metabolism. During aging, especially in barriques, exposure to oxygen and the lack of adequate topping up activate latent *Acetobacter bacteria*. These microorganisms are exceptionally resilient: they do not die during alcoholic fermentation and remain ready to convert ethanol into vinegar as soon as they encounter oxygen.

Cause of Formation	Typical Impact (g/L)	Operational Risk
Alcoholic Fermentation ( <i>S. cerevisiae</i> )	0.2 - 0.4 g/L	Low (Physiological process)
Malolactic Fermentation ( <i>O. oeni</i> )	0.4 - 1.0 g/L	Moderate (Monitor strains and times)
Deteriorated grapes / Botrytis	> 1.0 g/L	Critic
Management deficiencies (Oxygen)	Progressive	High (Degradation during refinement)

Table 1: Correlation between Formation Cause and Operational Risk

### 3. Comparative Analysis: OIV Official Method vs. CDR WineLab®

Analytical precision is the foundation of the technical validity of every winemaking decision. The choice of method determines not only the accuracy of the data, but also the efficiency of response times in the winery.

The **Official Method (OIV-MA-AS313-02)** It is a procedure based on steam distillation and subsequent titration with soda. Technically complex, it requires:

- After removal of CO<sub>2</sub> under vacuum.
- Mandatory correction for free and combined SO<sub>2</sub> by titration with **iodine**.
- Subtraction of sorbic acid and, in specific contexts, **salicylic acid** (used to stabilize samples), which requires an additional colorimetric protocol described in the OIV Annex.

The **CDR WineLab® System** offers an optimized solution for acetic acid determination: the application of highly sensitive LED photometric technology with intuitive methods using pre-filled and pre-calibrated reagents. Eliminating titration and sample handling eliminates the use of glassware and significantly reduces the use and exposure to hazardous solvents.

Criterion	Official OIV Method	CDR WineLab® System
<b>Procedural Complexity</b>	High (Distillation + multiple titrations)	Very Low (Few guided steps)
<b>Response Time</b>	30-45 minutes per sample	A few minutes
<b>Technical Accuracy</b>	Reference method (Titration)	Related to the reference (LED Photometry)
<b>Specialized Staff</b>	Essential (Chemical expertise)	Not required (Intuitive interface)
<b>Maintenance</b>	Constant cleaning and calibration	Absent
<b>Sustainability and Safety</b>	High production of chemical waste	Micro-quantity (Minimum impact/waste)

Table 2: Comparison of Analytical Methodologies

**Studies conducted by the University of Florence** confirm that **CDR WineLab® results are perfectly aligned with official methods**. This reliability, combined with automatic error checking, makes the system a professional precision instrument.

## CDR WineLab® System in the winery

Insourcing analyses eliminates the logistical barriers of external laboratories, eliminating waiting times and risks associated with sample transportation. This **analytical autonomy** allows for immediate corrective action, transforming real-time data into a competitive asset for the winery.

- **Decision-Making Autonomy and Timeliness:** Having data available within minutes allows for immediate interventions (cutting, racking, filtration) during harvest or maturation. This speed is essential for blocking bacterial growth detected at an early stage.
- **Simplicity and Operational Efficiency:** LED technology requires no maintenance or replacement like tungsten lamps. Thanks to the integrated "Help" interface, the system is accessible to all staff, ensuring repeatable results and reducing the risk of human error.
- **Sustainability and Safety:** The use of micro-quantities of reagents drastically reduces operator exposure to chemicals and minimizes waste, aligning with modern corporate sustainability standards.
- **Investment Versatility: In addition to wine, the system is configured to monitor vinegar production** with the same precision, offering unique operational flexibility for companies with diversified product lines.

## 5. Operational Summary for the Modern Producer

In a global market where winemaking precision is a prerequisite, a proactive approach to **acetic acid measurement** defines the boundary between controlled and managed production. Management success relies on the integration of three pillars: in-depth process knowledge, analytical precision, and technological efficiency. The integration of autonomous analysis systems like **CDR WineLab®** represents a strategic investment to transform quality control into a lasting competitive advantage. A simple, reliable, and rapid system can be a strategic partner for total control, from grape reception to bottling, elevating the excellence of your production.